

BACOBAN^{WB}

DISINFECTION & CLEANING CONCENTRATE

WATER-BASED

Alcohol-Free | Mixing table see reverse

Anti-Microbial effect for up to 10 days (Bacteria, Fungi and Viruses)

according with guideline 93/42/EEC; MDD

Composition:

100 g of solution containing: Benzalkonium chloride 26.0 g, sodium pyrithione 2.5 g, polycondensates, perfume substances, purified water, formaldehyde-free, phenol-free. Ingredients according to regulation (EC) No 648/2004: 15 % or over but less than 30 % cationic surfactants, disinfectants, perfumes (Linalool, Limone, Citral, Geraniol, Citronellol)

Efficacy and dosage	5 min	15 min	30 min	60 min	
VAH list (Bactericidal and yeasticidal efficacy)					
Tested in accordance with VAH requirements 2015 (with mechanical action, dirty conditions)	1,00 %				
EN 14476 and Guideline DVV / RKI					
Virucidal activity against enveloped viruses (incl. Coronavirus, HIV, HBV, HCV, Influenza, Herpes, Ebola virus...)	1,00 %				
EN 14476					
Adenovirus					1,00 %
Human papillomavirus (HPV)				1,00 %	
Norovirus			1,00 %		
Rotavirus	0,10 %				
Tested in accordance with EN guidelines					
EN 16615 (4-field-test)	1,00%				
EN 1040	0,25%				
EN 1275	0,25%				
EN 1276	0,7 %	0,50%			
EN 1650	0,50%	0,25%			
EN 13624 (C. albicans)	1,00%				
EN 13624 (T. interdigitale)	1,00%				
EN 13697			0,5 %		
EN 13727	1,00%				
EN 14561	3,00%			0,50%	
EN 14562				0,50%	

Mixing table	Final concentration					
	0,25 %	0,5 %	0,75 %	1 %	1,5 %	2 %
1 Litre	2,5 ml	5 ml	7,5 ml	10 ml	15 ml	20 ml
2 Litres	5 ml	10 ml	15 ml	20 ml	30 ml	40 ml
4 Litres	10 ml	20 ml	30 ml	40 ml	60 ml	80 ml
8 Litres	20 ml	40 ml	60 ml	80 ml	120 ml	160 ml
10 Litres	25 ml	50 ml	75 ml	100 ml	150 ml	200 ml

Danger



Residual antimicrobial effect for up to 10 days (tested according to ASTM E 2180).

Dermatologically tested (*at 1 %): excellent skin tolerance

Important!

The residual antimicrobial effect of Bacoban WB® is no excuse for not carrying out necessary disinfection measures in cases of visible and/or questionable contamination!

Use disinfectants safely. Always read the label and product information before use.

Safety instructions:

The Professional Association requires that suitable protective gloves be worn when handling disinfectants. Harmful if swallowed. Causes severe skin burns and eye damage. Very toxic to aquatic life. Keep out of reach of children. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. Dispose of contents/container to special waste collection point. Contains Linalool. May produce an allergic reaction. Toxic by eye contact. Pyridin-2-thiol-1-oxid, sodium salt. Use only with water and according to the mixing table. Do not use with other cleaners or disinfectants. Can be used on all surfaces. For possibly sensitive surfaces, we recommend testing for compatibility first on a less visible place. Store only in closed container, cool and dry. Container cannot be recycled until it is completely empty. Do not use product after expiry date. Dispose of properly in accordance with local and official regulations.

The ready-to-use solution based on the mixing table does not have to be labelled with a danger symbol.



BACOBAN^{DL} 1% DISINFECTION & CLEANING READY-TO-USE

WATER-BASED

Alcohol-Free | For surfaces, medical fittings and equipment

Anti-Microbial effect for up to 10 days (Bacteria, Fungi and Viruses)

according with guideline 93/42/EEC; MDD

Ingredients:

100 g of solution containing: Benzalkonium chloride 0,26 g, sodium pyrithione 0,025 g, polycondensates, perfume substances, purified water. Ingredients according to regulation (EC) No 648/2004: less than 5 % cationic surfactants, disinfectants, perfumes.

Intended use:

Ready-to-use disinfectant product with residual effect. Does not contain Aldehyde or Phenol.

Efficacy and application:

Tested in accordance with EN guidelines:

EN 1040: 5 min
 EN 1275: 5 min
 EN 1276: 5 min
 EN 1650: 5 min
 EN 13697: 15 min
 EN 13727: 5 min
 EN 13624: 5 min (C. albicans)
 EN 13624: 5 min (T. interdigitale)
 EN 14561: 60 min
 EN 14562: 60 min
 EN 16615: 5 min (4-field-test)
 Virucidal activity against enveloped viruses according to EN 14476 under dirty conditions: 5 min (incl. Coronavirus, HIV, HBV, HCV, Influenza, Herpes, Ebola virus,...)
 EN 14476: Rotavirus 5 min
 Norovirus 30 min
 Human papillomavirus 60 min

Residual antimicrobial effect for up to 10 days (tested according to ASTM E 2180).

Dermatologically tested (epicutaneous test with 30 test persons): excellent skin tolerance



Safety instructions:

Wear gloves when using disinfectant products. Keep out of reach of children. Can be used on all surfaces. For possibly sensitive surfaces, we recommend testing for compatibility first on a less visible place. To be kept in a dry, aired place. Close properly after use. Recycle and dispose of the container in accordance with current legislation, preferably via an approved collector or company. Do not use after the sell-by date. A safety sheet can be obtained on request.

Important!

The residual effect of Bacoban® DL is no excuse for not carrying out necessary disinfection measures in cases of visible and/or questionable contamination!

Use disinfectants safely. Always read the label and product information before use.

Directions for use: Spread Bacoban® DL evenly over the surface with a cloth and leave to dry. Do not rinse.

BACOBAN^{DL} 3% DISINFECTION & CLEANING READY-TO-USE

WATER-BASED

Alcohol-Free | For surfaces, medical fittings and equipment

Anti-Microbial effect for up to 10 days (Bacteria, Fungi and Viruses)

according with guideline 93/42/EEC; MDD

Composition:

100 g solution contains: Benzalkonium chloride 0.78 g, sodium pyrithione 0.075 g, polycondensates, fragrances, purified water

Intended purpose:

Ready-to-use disinfection solution for use in connection with an Adexano nebuliser with long-term antimicrobial effect. Formaldehyde-free, phenol-free

Efficacy and application:

Tested acc. the DGHM (German Society for Hygiene and Microbiology) guideline (high organic load) 1.0 %, 5 min.; 0.75 %, 15 min.
Tested acc. EN guidelines:

DIN EN 1040: 5 min.
DIN EN 1275: 5 min.
DIN EN 1276: 5 min.
DIN EN 1650: 5 min.
DIN EN 13697: 5 min.
DIN EN 14476: 2005 Adenovirus 2% 240 min.
DIN EN 14476: 2005 Rotavirus 0.1% 5 min.
DIN EN 14476: 2007-2 Norovirus 2% 240 min.



Safety guidelines:

Do not stay in the same room whilst disinfection is being carried out with the Adexano nebuliser. Keep windows and doors closed. Do not swallow or inhale. Avoid eye contact. Rinse affected areas thoroughly with water after contact with eyes or mucous membranes and contact your doctor. Keep away from children.

Can be used on all surfaces. We recommend that you test sensitive surfaces in an inconspicuous place with Bacoban®DL 3 %. The product should be kept only in a closed vessel. The container may only be submitted for recycling in a completely empty condition. Do not use after the expiry date. Dispose of in a suitable manner, taking the local and official specifications into account

Important! The long-term antimicrobial effect of Bacoban® DL 3 % does not mean that necessary disinfection measures for visible and/or questionable contamination should be neglected.

Use disinfectants safely. Always read the marking and the product information prior to use.